

Package ‘Affyhgu133A2Expr’

April 13, 2017

Type Package

Title Affymetrix Human Genome U133A 2.0 Array (GPL571) Expression Data Package

Version 1.10.0

Date 2014-2-13

Author Zhicheng Ji, Hongkai Ji

Maintainer Zhicheng Ji <zji4@jhu.edu>

Description Contains pre-built human (GPL571) databases of gene expression profiles. The gene expression data was downloaded from NCBI GEO and preprocessed and normalized consistently. The biological context of each sample was recorded and manually verified based on the sample description in GEO.

License GPL (>=2)

Depends R (>= 2.10)

biocViews Genome, Homo_sapiens_Data, GEO

NeedsCompilation no

R topics documented:

Affyhgu133A2Expr-package	1
Affyhgu133A2Expr	2
Affyhgu133A2Exprtab	3

Index	5
--------------	----------

Affyhgu133A2Expr-package
Affymetrix Human Genome U133A 2.0 Array (GPL571) Expression Data Package

Description

Contains gene expression profiles from Affymetrix Human Genome U133A 2.0 Array (GPL571). This package is primarily designed for GSCA (Gene Set Context Analysis). All gene expression data are downloaded from NCBI GEO. Gene expression data were preprocessed and normalized consistently using fRMA. Gene expression values are further standardized using gene expression barcode approach. The biological context of each sample was recorded and manually verified based on the sample description in GEO. Gene expression profiles are stored as hdf5 format.

Details

Package: Affyhgu133A2Expr
Type: Package
Version: 1.0.0
Date: 2014-4-9
License: GPL 2.0

Author(s)

Author: Zhicheng Ji, Hongkai Ji Maintainer: Zhicheng Ji <zji4@jhu.edu>

References

- McCall M.N., Bolstad B.M., and Irizarry R.A. (2010) Frozen robust multiarray analysis (fRMA). *Biostatistics* 11, 242-253.
- McCall, M. N., Uppal, K., Jaffee, H. A., Zilliox, M. J., & Irizarry, R. A. (2011). The Gene Expression Barcode: leveraging public data repositories to begin cataloging the human and murine transcriptomes. *Nucleic acids research*, 39(suppl 1), D1011-D1015.
- Barrett T., et al. (2007) NCBI GEO: mining tens of millions of expression profiles - database and tools update. *Nucl. Acids Res.* 35, D760-D765.

Affyhgu133A2Expr *Data of human gene expression profiles from the Affymetrix Human Genome U133A 2.0 Array (GPL571).*

Description

The data set contains 313 human profiles on 12494 genes downloaded from NCBI GEO. Gene expression data were preprocessed and normalized consistently using fRMA. Gene expression values are further standardized using gene expression barcode approach. The biological context of each sample was recorded and manually verified based on the sample description in GEO. The gene expression value matrix is stored in hdf5 format using rhdf5 package.

Details

This data package contains expression values of 12494 genes and 313 samples measurements from NCBI GEO obtained using the GPL571 platform. Gene expression data were preprocessed and normalized consistently using fRMA. Gene expression values were further standardized using gene expression barcode approach. Probeset with the largest coefficient of variation from all probesets that corresponds to the same gene is retained, so that each gene uniquely matches to one row in the database. The biological context of each sample was also recorded and manually verified based on the sample description in GEO. To enhance the reading speed of the dataset, the gene expression value matrix is stored in a hdf5 format using rhdf5 package. The rows of the matrix represents samples and the columns of the matrix represent genes. Notice that all values are 1000 times the actual value so that the values can be stored as integers to minimize file size and reading time. The package is specifically designed to be manipulated by GSCA package so users are not expected to read the expression values by themselves. The sample id, sample type, and experiment id for each sample in the gene expression compendium are also included in this data package.

Source

www.ncbi.nlm.nih.gov/geo/

References

McCall M.N., Bolstad B.M., and Irizarry R.A. (2010) Frozen robust multiarray analysis (fRMA). *Biostatistics* 11, 242-253.

McCall, M. N., Uppal, K., Jaffee, H. A., Zilliox, M. J., & Irizarry, R. A. (2011). The Gene Expression Barcode: leveraging public data repositories to begin cataloging the human and murine transcriptomes. *Nucleic acids research*, 39(suppl 1), D1011-D1015.

Barrett T., et al. (2007) NCBI GEO: mining tens of millions of expression profiles - database and tools update. *Nucl. Acids Res.* 35, D760-D765.

Affyhgu133A2Exprtab *Reference table for Affyhgu133A2Expr gene expression compendium*

Description

Contains the sample id, sample type, and experiment id for each sample in the Affymetrix Human Genome U133A 2.0 Array (GPL571) gene expression compendium.

Usage

```
data(Affyhgu133A2Exprtab)
```

Format

A data frame with 313 observations on the following 3 variables.

SampleID a character vector

ExperimentID a character vector

SampleType a character vector

Details

SampleID is a GSM ID that NCBI GEO uses to as a sample identifier. ExperimentID is a GEO ID that NCBI GEO uses to identify an experiment. SampleType denotes the cell type or tissue and whether the sample is given a specific treatment or in a specific condition.

Source

www.ncbi.nlm.nih.gov/geo/

References

McCall M.N., Bolstad B.M., and Irizarry R.A. (2010) Frozen robust multiarray analysis (fRMA). *Biostatistics* 11, 242-253.

McCall, M. N., Uppal, K., Jaffee, H. A., Zilliox, M. J., & Irizarry, R. A. (2011). The Gene Expression Barcode: leveraging public data repositories to begin cataloging the human and murine transcriptomes. *Nucleic acids research*, 39(suppl 1), D1011-D1015.

Barrett T., et al. (2007) NCBI GEO: mining tens of millions of expression profiles - database and tools update. *Nucl. Acids Res.* 35, D760-D765.

Examples

```
## Load the reference table  
data(Affyghu133A2Exprtab)  
str(Affyghu133A2Exprtab)
```

Index

*Topic **datasets,**
 Affyhu133A2Exprtab
 Affyhu133A2Exprtab, [3](#)

*Topic
 datasets,GPL571,database,Affyhu133A2Expr
 Affyhu133A2Expr, [2](#)

*Topic **package, database,**
 Affyhu133A2Expr
 Affyhu133A2Expr-package, [1](#)

Affyhu133A2Expr, [2](#)
Affyhu133A2Expr-package, [1](#)
Affyhu133A2Exprtab, [3](#)